



NEUROLINE
PRODUCT CATALOGUE
2023

We reserve the right to change specifications, features and solutions without notice.

PLEASE NOTE THAT SOME ITEMS IN THIS CATALOGUE MAY HAVE BEEN DISCONTINUED OR TAKEN OUT OF PRODUCTION. PRICES, COLOURS AND/OR COMPOSITION OF OUR PRODUCTS

MAY CHANGE OVER TIME. WE RECOMMEND A VISIT TO OUR WEBSITE FOR UPDATED VERSIONS OF OUR BROCHURES.

© 2023 EMD KFT. ALL RIGHTS RESERVED. NO PART OF THIS PUBLICATION MAY BE REPRODUCED, STORED IN A RETRIEVAL SYSTEM, OR TRANSMITTED IN ANY FORM OR BY ANY MEANS,

ELECTRONIC, MECHANICAL, PHOTOCOPYING, RECORDING OR OTHERWISE, WITHOUT PRIOR WRITTEN CONSENT OF EMD KFT.

About Us

EMD PROFILE

COMPANY HISTORY

Neurosurgery
NEUROSMART HIGH-SPEED MOTOR SYSTEM

20 Craniotomy
HANDPIECES
CUTTERS

26 Attachments

STANDARD HANDPIECES
ANGLED HANDPIECES
SLIM ATTACHMENTS
BALLTIP BURRS

TRANSFORAMINAL

Endoscopic Spine Surgery

HANDPIECES OF ENDOSCOPES

STENOSIS
INERLAMINAL I. - II.

48 Cranial Perforation

PERFORATOR HANDPIECE
DISPOSABLE CRANIAL PERFORATORS

REUSABLE CRANIAL PERFORATORS

Accessories

ACCESSORIES OF IRRIGATION
WASH-RINSE DEVICE
STORAGE CONTAINER
SERVICE OILING SPRAY



Be Ahead Of The Curve With Our Latest Innovations

WE ARE FAMILY OF MOTIVATED INDIVIDUALS, DRIVEN BY
OUR DESIRE TO MAKE INNOVATIVE, EXCELLENT QUALITY
MEDICAL PRODUCTS AVAILABLE TO MODERN MEDICINE IN
A SUSTAINABLE WAY.

At EMD we offer quality design, production and services, that enable medical professionals to perform advanced minimal invasive operations.

Our carefully accumulated knowledge and expertise enable us to create singular tools and devices that medical professionals can use sustainably. To offer our customers the best of the best, our products are manufactured with meticulous care from high quality materials.

We are dedicated and passionate in everything we do. Through our performance and solution driven mindset, we have the opportunity to make a difference. We are innovation and business oriented with a strong belief in daring to do things better.

From CEO

66

Everything we do is rooted in the belief that, together with our colleagues, teams and partners, we can create positive change to support modern medicine through our high-end products and services.

Each time before we embark on a new development, we have to make sure we possess the essential attributes that will fuel us to work for years to come for the very moment we first hold a new product in hand. These are: know-how, perseverance, proof and faith.

We've been working hard and have come a long way over the last 25 years. When EMD was founded in 1997, we were thinking about how we could design and manufacture minimal invasive medical devices. Since then, we've made such a progress as a team that today, however, we are thinking about how to become one of the best and most innovative companies in the industry.

We are proud of the progress we're making, and are ensuring that our expertise and competencies are essential elements of EMD's future that we will continue to embrace and grow."





8 Product Catalogue

About Us **emd** 9

EMD is founded as a repair service provider. Main activity is repair of endoscopic optics and endoscopic 1997 instruments. We launch our laparoscopic system being first own product 2001 development. EMD steps out of being a local supplier and starts export 2005 activities in various markets. EMD begins to develop and produce neurosurgery devices. EMD relocates to new headquarters with more on site facilities. 2012 Company profile continues to expand. The company obtains CE approval for own neurosurgical and 2017 laparoscopy development. Product portfolio is further increased by six additional EMD 2021 product group receiving CE mark. First product group got FDA approval. 2022

Our Vision

We believe in leading not following. We turn fault into our virtue and trust our nevercompromising nature dictating our rhythm, which takes us to innovative ideas and exceptional product developments.

We are convinced that possibility is bigger than the fear of failure. That is why we dare to dream big and are not afraid to be challenged on the way to success. Our vision is to be industry leaders in medical device development and manufacture.

Our Mission

OUR MISSION IS TO SUPPORT MODERN MEDICINE BY OFFERING ADVANCED, FIRST-CLASS QUALITY MEDICAL DEVICES.

Our Values

We like to walk our talk and lead by example. Our values are part of who we are, what we stand for and how we act.

Our great colleagues make the difference. It's when we share our skills, knowledge and expertise, we become one team.

We encourage and help each other to achieve our goals, while always having our customers' and their patients' best interest in mind.

Markets & Expansion

EMD is expanding, with a focus on multi-channel sales. We treat and nurture long term partnerships and business growth as one of our top priorities. We offer advanced and unique medical tools for minimal invasive surgeries that want to give professionals unbeatable value with the best combination of innovation, quality, price and sustainability.

When it comes to expansion strategies and activities, we spare no time to find the most prosperable opportunities we are able to find out there, that not only benefit us, but our partners too. With the continuous product developments we have been making recently, we are enthusiastic to be on the path greatly expanding our existing presence on the market.

Research & Development

Our network of medical professionals, together with the global industry needs and challenges set new directions to our developments.

Through our product developments, our priority is to effectively support surgical procedures, which are still very difficult or can only be performed with conventional surgery, so that they could be carried out with as few and little incisions as possible. In the near future, we continue to focus on minimal invasive developments and present new products for spine,neuro,- and laparoscopic surgery.

Innovative Solutions

Innovative product developments enable us to find and implement new technologic solutions, which are not yet in practice.

At the birth of a new idea, we start exploring ways of execution and won't stop until we find the right technology that can be applied for best results. This way, we are inspired to never stand still, we find new discoveries and are able to create that was thought impossible.

Professional Experience

COMBINED WITH BRILLIANT IDEAS IN THE SERVICE OF NEUROSCIENCE.

NEUROSURGERY IS ABOUT PRECISION AND CAREFUL PLANNING.

WITH THE BEST OF OUR KNOWLEDGE AND EXPERIENCE, IN CLOSE CO-OPERATION WITH MOST RENOWNED NEUROSURGEONS, WE OFFER OWN DEVELOPED, PATENTED PRODUCTS IN THE FIELDS OF CRANIAL AND SPINE SURGERY.



With **MULT:SMART**

WE OFFER MULTI-CHOICE SOLUTIONS

DESIGNED FOR ALL NEUROSURGICAL APPLICATIONS, NEUROSMART IS AN ADVANCED ELECTRIC MOTOR SYSTEM OF EMD'S NEUROLINE PRODUCT LINE WITH ITS EXTENDED USABILITY FROM CONVENTIONAL SKULL AND BONE SURGERY TO MINIMAL INVASIVE OPERATIONS.

In addition to the High Speed Mode which is most commonly used for skull and bone surgery, the MULTISMART-featured control unit has been enhanced with two additional functions (Endoscopic Spine & Slim Attachment Modes) that facilitate the system's applicability to spinal surgery and thinner, lower-speed handpieces.

Intelligent improvements eg. Automatical Handpiece Detection and Speed Range Optimalisation of the Dedicated Modes maximize patient safety and support user experience.

The upgraded control unit is compatible with the NSS-120 HS engine, which already delivers



High Speed RPM MAX 80 000



Endoscopic Spine RPM MAX 20 000



Slim Attachment RPM MAX 60 000



Maximum Safety. High Power.

- MULTI-APPLICATION MODES OF HIGH SPEED/ENDOSCOPIC SPINE/SLIM ATTACHMENT
- **OPTIMAL SPEED RANGES GRANT PATIENT SAFETY**
- **OUTSTANDING MOTOR PERFORMANCE OF 120W**
- SPEED RANGE 5.000-80.000 RPM
- COST-EFFECTIVE SOLUTION ONE MOTOR OF ALL ATTACHMENTS
- AUTOMATIC HANDPIECE DETECTION FOR MAXIMUM SAFETY IN HIGH SPEED MODE
- INTEGRATED IRRIGATION PUMP WITH 5-DEGREE ADJUSTMENT
- TOUCH SCREEN DISPLAY FOR EASY HANDLING AND CLEANING
- ONE-CLICK CONNECTION OF SYSTEM ELEMENTS



With its one-piece glass front panel, the NeuroSmart Control Unit was designed with medical hygiene and user practicality in mind.

The touch screen display offers easy handling and all-inclusive cleaning options. The user-friendly menu supports easy navigation in functions and settings.





LIGHT/DARK **APPEARANCE**

The user has the possibility to switch between light and dark view modes at any point during procedure.

INTEGRATED IRRIGATION PUMP

During neurosurgical procedures, it is essential to ensure about proper cooling at surgical area. The Control Unit features integrated irrigation pump with variable flow control.







MULT SMART High Speed

The intelligent, **emd**®-engineered Motor System recognizes whether a Perforator Handpiece or high speed Handpiece is connected to the motor.

At the same time as the Perforator Handpiece is connected by the user, the Control Unit automatically sets for its most optimal speed for **emd** Cranial Perforators, preventing inappropriate speed starting.

For maximum patient protection, speed drive is automatically optimized to ensure perforation starts at appropriate rpm. Therefore, with the High-Speed Motor System the risk of an accidental high speed start is completely avoided. When the adapter is disconnected, the control unit resets to latest speed use of High Speed Mode.

AUTOMATIZED MODE SETUP

Connection and disconnection of the Perforator Handpiece results in an immediate, automatized change of user modes and display functions.

In Perforation Mode speed is fixed on 1.000 RPM.

Irrigation can be activated/deactivated by using the rinse button of the footswitch expect for use in Perforation Mode.





MULT:SMART Endoscopic Spine

For minimal invasive applications in Neurosurgery, dedicated handpieces are developed for the NeuroSmart Motor System.

In Endoscopic Spine Mode the controller software sets the maximum speed that can be adjusted to the limit specified for the Spine Handpieces. This feature eliminates incorrect, excessive speed settings, and provides safety for both user and patient, while protecting the handpiece from premature failure.

Speed is limited in range of 3.000-20.000 RPM. Irrigation is disabled in this mode.



MULT:SMART Slim Attachment

The purpose of Slim Attachment Mode is to ensure safe, error-free use for the user when applying slim attachments with the NeuroSmart Motor System This upgrade mode of operation is designed to allow the use of slim handpieces with limited speed range.

Speed is limited in range of 10.000-60.000 RPM. Irrigation is disabled in this mode.





SMART & SUPERB POWERED MOTOR

Highly cost-efficient, due to its high-excursion motor power, the NeuroSmart High-Speed Motor by itself manages to drive the entire range of **emd** Handpiece attachments at appropriate speed. Using One-Click Technology, the complete range of Hand,- and Nosepiece tools are easily connected to the motor by a single snapping movement.

Integrated bend relief to prevent cable torsion and increase protection of connection.



office

CABLE CONNECTED FOOTSWITCH

The combined Footswitch pedal is used to start/stop motor operation, while using the push button, rinsing mode is turned on/off.

The footswitch construction complies with IPX8 standard and is therefore easy to clean and disinfect.

Pedal and push button positioning delivers comfort even during longer user times.

Standard Set Configuration



NSS-110

CONTROL UNIT

Power consumption: 300 VA Power Supply: 110-240 VAC; 50-60 Hz Irrigation capacity: 0-50 ml/min. Weight: 5kg Device class: Ilb Protection class: I



NSS-120

HIGH SPEED MOTOR with integrated cable 3m long connecting cable

Autoclavable Max speed: 80.000 rpm Output power: 120W Device class: IIb



NSS-130

DOUBLE FOOTSWITCH 4m long connecting cable

Device class: Ilb Standard IPX8



20 Product Catalogue 21

WE BELIEVE THAT

Functionality And Design

CREATE WORK OF ART.

EITHER WE ARE WORKING ON A NEW PRODUCT DEVELOPMENT OR WE ARE RETHINKING AN OLDER DESIGN, OUR TOP PRIORITY IS CREATING HIGH-END PRODUCTS THAT COMPLY WITH THE CHALLENGES OF THE TIME, BOTH IN TERMS OF FUNCTION AND USABILITY.



Absolute Performance. Without Compromises.

EACH **emd** HANDPIECE COLLECTION IS MADE OF STAINLESS STEEL, GUARANTEEING LONG PRODUCT LIFE. WITH THE SLEEKEST DESIGN ON THE MARKET, A CLEAN SHAPE AND SURFACE MAKES OUR CRANIOTOME HANDPIECE ESPECIALLY EASY TO CLEAN AND HOLD. THE ERGONOMIC PEN-GRIP SOLUTION ALLOWS COMFORTABLE USE AND A STABLE GRIP DURING OPERATION.

The One-Click Technology provides easy connection of the Craniotome Handpieces by simply snapping it to the motor.

Both the inserted craniotome cutter and duraguard are simultaneously locked in position with rotary knob. For greatest safety, starting of the motor is not possible until rotary knob is locked.

For more user support, duraguards and craniotome cutters are marked with size identifiers. Duraguards are available in fixed and rotatable configurations.

The range also includes additional models for special application areas. Eg. special head for craniotomy nearby the sinus or duraguards with increased head length for laminotomy.

Craniotome/laminotome heads with fix duraguard are also available with irrigation. The feature ensures a consistent low temperature during cutting.









Illustration of emd[®]Craniotome Handpiece



Illustration of **emd**[®]Craniotome Handpiece with Craniotome Head



NSH-320

NSH-332

CRANIOTOME HEAD Fix - Without Duraguard

NLC-1T 3.5mm exposure NLC-1ST 3.5mm exposure NLC-2T 5.5mm exposure NLC-2ST 5.5mm exposure



NSH-331R

CRANIOTOME HEAD, L11 R Rotatable - Duraguard - Size 1 NLC-1T NLC-1ST



NSH-332R

CRANIOTOME HEAD, L15 R Rotatable - Duraguard - Size 2 NLC-2T NLC-2ST



NSH-333R

CRANIOTOME HEAD, L25 R Rotatable- Duraguard - Size 3

NLC-3T NLC-3ST



NSH-334R

LAMINOTOME HEAD, L35 R Rotatable - Protection Guard (max. speed 30,000 rpm) Use with cutters

NLC-4T NLC-4ST



NSH-331

CRANIOTOME HEAD, L11 Fix - Duraguard - Size 1

NI C-1T NI C-1ST



NSH-332

CRANIOTOME HEAD, L15 Fix - Duraguard - Size 2 Use with cutters:

NLC-2T NLC-2ST



NSH-333

CRANIOTOME HEAD, L25 Fix - Duraguard - Size 3 NLC-3T

NLC-3ST



NSH-334

LAMINOTOME HEAD, L35 Fix - Protection Guard (max. speed 30,000 rpm) Use with cutters NLC-4T NLC-4ST



SPEED RANGE 10.000-80.000 RPM (EXCEPT FOR LAMINOTOME)

emd CRANIOTOME HANDPIECES ARE DESIGNED FOR BOTH ADULT AND PEDIATRIC SURGERY

emd OFFERS THREE FOOTED DURAGUARDS; ONE LAMINOTOME HEAD; ONE SPECIALTY HEAD FOR SINUS

QUICK AND EASY RELEASE CHUCK OF DURAGUARDS

CRANIOTOME/LAMINOTOME HEADS WITH FIX DURAGUARD ARE ALSO AVAILABLE WITH IRRIGATION

PRECISELY DESIGNED EDGE CONSTRUCTIONS, DURABLE BLADES ARE THE MAIN FEATURES OF NEUROLINE CRANIOTOME CUTTERS.

AVAILABLE IN FOUR STANDARD LENGTHS TO FIT THE RIGHT SIZE DURAGUARD OR EVEN A SPECIAL CRANIOTOME HEAD.

THREE CUTTING BLADES DELIVER THE BEST PERFORMANCE BOTH FOR STRAIGHT AND SPIRALED MODELS.



NLC-1T

CRANIOTOME CUTTER, L11 Straight with 3 cutting edge - Size 1 Use with duraguard: NSH-331 NSH-331R

NLC-2T

CRANIOTOME CUTTER, L15 Straight with 3 cutting edge - Size 2 Use with duraguard: NSH-332 NSH-332R

NLC-3T

CRANIOTOME CUTTER, L25 Straight with 3 cutting edge - Size 3 NSH-333R

NLC-4T

CRANIOTOME CUTTER, L35 Straight with 3 cutting edge - Size 4 NSH-334 NSH-334R



NLC-1ST

CRANIOTOME CUTTER, L11 Spiral with 3 cutting edge - Size 1 Use with duraguard NSH-331 NSH-331R

NLC-2ST

CRANIOTOME CUTTER, L15 Spiral with 3 cutting edge - Size 2 Use with duraguard: NSH-332R

NLC-3ST

CRANIOTOME CUTTER, L25 Spiral with 3 cutting edge - Size 3 Use with duraguard: NSH-333 NSH-333R

NLC-4ST

CRANIOTOME CUTTER, L35 Spiral with 3 cutting edge - Size 4 Use with duraguard: NSH-334 NSH-334R



26 Product Catalogue 27

WE BELIEVE THAT

Uncompromising Quality

IS A MATTER OF HIGH EXPECTATIONS.

GUARANTEEING THE RELIABILITY OF AN **emd**®PRODUCT REQUIRES THE APPLICATION OF A MULTITUDE OF SKILLS. FROM THE PRECISE SELECTION OF RAW MATERIALS AND MACHINING TECHNOLOGIES TO THE COMPLEXITY OF PRODUCT CERTIFICATION, EVERY STEP IS TAKEN TO ENSURE THAT OUR QUALITY ITEMS MEET THE HIGHEST SAFETY STANDARDS FOR BOTH THE PATIENT AND THE SURGEON.



Absolute Performance. Without Compromises.

EACH **emd** HANDPIECE COLLECTION IS MADE OF STAINLESS STEEL, GUARANTEEING A LONG PRODUCT LIFE. FEATURING A PURPOSEFULLY SLEEK DESIGN, THE CLEAN SHAPE AND SURFACE MAKES OUR STANDARD HANDPIECES ESPECIALLY EASY TO CLEAN. THE ERGONOMIC PEN-GRIP SOLUTION ALLOWS COMFORTABLE USE AND A STABLE GRIP DURING OPERATION.

The One-Click Technology provides easy connection of the Standard Handpieces by simply snapping it to the motor.

Both the inserted tool cutter and the handpiece are simultaneously locked in position with rotary knob. For greatest safety, starting of the motor is not possible until rotary knob is locked.

The special Tool Grasping Technology of attachments supports product selection and economic operation.







Features •

\$\ightharpoonup \text{STANDARD HANDPIECES ARE COMPATIBLE WITH BURRS \(\text{\rho} \) 2,3MM

STANDARD HANDPIECES OF THE NEUROLINE FAMILY

UNIFORM-SIZED ADJUSTABLE EXPOSURE TOOLS ARE COMPATIBLE WITH ALL



Burrs With Benefits

USER-FRIENDLY AND COST-EFFICIENT HS TOOLS ARE AVAILABLE WITH THE MOST COMMONLY USED GEOMETRIES.

THANKS TO THE SPECIAL TOOL GRASPING TECHNOLOGY, ALL STANDARD AND ANGLED **emd**® NEUROLINE HANDPIECES ARE TO BE USED WITH UNIFORM-SIZED BURRS AND TOOLS.




ONE SIZE. MANY APPLICATIONS.



EXPOSURE IN 8 POSITIONS

Desired tool length is easily adjustable according to user preferences.



Absolute Performance. Without Compromises.

EACH **emd** HANDPIECE COLLECTION IS MADE OF STAINLESS STEEL, GUARANTEEING A LONG PRODUCT LIFE. FEATURING A PURPOSEFULLY SLEEK DESIGN, THE CLEAN SHAPE AND SURFACE MAKES OUR SLIM HANDPIECES ESPECIALLY EASY TO CLEAN. THE ERGONOMIC PEN-GRIP SOLUTION ALLOWS COMFORTABLE USE AND A STABLE GRIP DURING OPERATION.

Several neurosurgical applications require more precise access to the operating area.

Slim attachments have a detachable construction of angeled handpiece and a small bore nosepiece. The surgeon can simply exchange the slim head on the handpiece to the suitable length or form model without any tools needed. Connection between the head and the handpiece is an **emd**[®] patent.

Slim Attachments are mainly used in spine surgery applications. All handpiece and nosepiece combinations are well balanced which is a great benefit during longer operations.



MULT:SMART





NSH-510

SLIM HANDPIECE Angled



NSH-521

LINEAR SLIM HEAD Working length 85mm, ø 4.7mm



NSH-522

LINEAR SLIM HEAD Working length 105mm, ø 4.7mm



NSH-523

LINEAR SLIM HEAD Working length 125mm, ø 4.7mm



NSH-531

CURVED SLIM HEAD Working length 85mm, ø 4.7mm



NSH-532

CURVED SLIM HEAD Working length 105mm, ø 4.7mm



NSH-533

CURVED SLIM HEAD Working length 125mm, ø 4.7mm





SLIM HEADS ARE COMPATIBLE WITH BURRS Ø 1.0MM

EACH SLIM ATTACHMENT CAN BE ROTATED AND LOCKED IN THE DESIRED POSITION IN A 360° RANGE












\mathbb{Z}	2	185	NLB-Z-A110A
(MM = Ø)	3		NLB-Z-A120A
	4		NLB-Z-A130A
\mathbb{Z}	2	105	NLB-Z-A110B
	2	_	NLB-Z-A120B
٥	4		NLB-Z-A130B
\mathbb{Z}	2	1125	NLB-Z-A110C
(MM = Ø)	3		NLB-Z-A120C
٣	Δ		NLB-Z-A130C



General	
eatures	

SPEED RANGE 10.000-60.000 RPM

BURRS SHOWED ON THIS PAGE ARE COMPATIBLE WITH **emd** NSH-5XX SLIM HEADS

SINGLE USE PRODUCT SUPPLIED IN STERILE PACKING



N N	2	L85	NLB-Z-A140A
= 0	3		NLB-Z-A150A
٥	4		NLB-Z-A160A
(N	2	105	NLB-Z-A140B
(MM = Ø)	3	_	NLB-Z-A150B
ی	4		NLB-Z-A160B
(N	2	125	NLB-Z-A140C
(Ø = MM)	3	_	NLB-Z-A150C
S.	4		NLB-Z-A160C



SPEED RANGE 10.000-60.000 RPM

BURRS SHOWED ON THIS PAGE ARE COMPATIBLE WITH **emd** NSH-5XX SLIM HEADS

SINGLE USE PRODUCT SUPPLIED IN STERILE PACKING

36 Product Catalogue 37

WE SHOW MORE OF OUR Expertise in Endoscopy

MODERN SURGICAL TECHNOLOGIES ARE GAINING GROUND IN AN INCREASING

NUMBER OF APPLICATIONS. AT **emd**® WE ARE COMMITTED TO MAKING OUR

DECADES OF DEVELOPMENT AND MANUFACTURING EXPERIENCE IN ENDOSCOPY

AVAILABLE TO AS MANY SURGICAL FIELDS AS POSSIBLE.



Surger

ENDOSCOPIC SPINE

ENDOSCOPY IS AN ADVANCED TECHNIQUE IN OPERATION PROCEDURES, WHICH BRINGS THE ADVANTAGES OF THE MINIMAL INVASIVE SURGERY METHOD IN MANY MEDICAL FIELDS.

SPINE SURGERY IS A TYPICAL APPLICATION OF MIS IN NEUROSURGERY THAT CAN BENEFIT THE MOST FROM THE ADVANTAGES OF THIS SURGICAL METHOD, EITHER FROM THE PATIENT'S OR THE PHYSICIAN'S POINT OF VIEW.

For endoscopic spine surgery, disposable handpieces with and without protection and with distal and beveled protection guide tubes have been developed for this purpose. Single use spine handpieces have a compact design containing the burrs available with four burr shapes in two working lengths. With this economic construction spine handpieces can be connected directly to the motor ignoring the demand for any additional attachment of the motor.

emd Endoscopes and Handpieces designed to perform endoscopic spinal surgery using the Stenosis, Interlaminal and Transforaminal approaches are designed with some overlap in size, allowing the user to minimize the cost of the operative devices.





NET-100A15

ENDOSCOPE

Working channel: 7.1mm Irrigation-suction channel: 2 x 1.7mm Working length: 125mm Angle of view: 15°

NED-100AA

DILATOR, SINGLE HOLE

Length: 170mm

NES-100AA

WORKING TUBE, BEVELED

nner diameter: ø 10.1mm Length: 126mm

NES-100AB

S

ENO

WORKING TUBE, FLAT

Outer diameter: ø 12mm Inner diameter: ø 10.1mm Length: 126mm

NET-070A30

ENDOSCOPE

Outer diameter: ø 6.9mm Working channel: 4.2mm Irrigation-suction channel: 2 x 1.5mm Working length: 130mm Total length: 204mm Angle of view: 30°

NED-070AB

DILATOR, 2 CHANNEL

Length: 175mm

NES-070AA

WORKING TUBE, BEVELED

Outer diameter: ø 9mm Inner diameter: ø 7mm Length: 131mm

ENDOSCOPE

NET-070B30

Outer diameter: ø 6.9mm Working channel: 4.1mm Irrigation-suction channel: 2 x 1.4mm Working length: 181mm Total length: 255mm Angle of view: 30°

NED-070BB

DILATOR, SINGLE HOLE

Length: 226mm

AMINA

SFOR/

NES-070BA

WORKING TUBE, BEVELED

nner diameter: ø 7mm Length: 182mm

NES-070BC

WORKING TUBE, ELEVATED

Inner diameter: ø 7mm Length: 182mm



AL

ERLAMIN

NET-070B30 **ENDOSCOPE**

Outer diameter: ø 6.9mm Working channel: 4.2mm Irrigation-suction channel: 2 x 1.4mm Working length: 181mm Total length: 255mm Angle of view: 30°

NED-070BB

DILATOR, 2 CHANNEL

Diameter: ø 7mm Length: 226mm

NES-070BA

WORKING TUBE, BEVELED

Inner diameter: ø 7mm Length: 182mm



NLB-Z-A210A

SPINE HANDPIECE, L220, DIAMOND **BALLTIP BURR IN SUPPORT TUBE**

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A220A

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 3mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A230A

SPINE HANDPIECE, L220, DIAMOND

BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mn Ball diameter: ø 4mm

Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A240A

SPINE HANDPIECE, L300, DIAMOND

BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 4mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A250A

SPINE HANDPIECE, L220, DIAMOND

BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 5mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A260A

SPINE HANDPIECE, L300, DIAMOND **BALLTIP BURR IN SUPPORT TUBE**

Ball diameter: ø 5mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A270A

SPINE HANDPIECE, L220, DIAMOND

BALLTIP BURR IN SUPPORT TUBE

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A280A

SPINE HANDPIECE, L300, DIAMOND

BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 6mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A010A

SPINE HANDPIECE, L220, STAINLESS STEEL

BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 3mm

Working length: 220mn

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A020A

SPINE HANDPIECE, L300, STAINLESS STEEL

BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 3mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A030A

SPINE HANDPIECE, L220, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mr

Ball diameter: ø 4mm

Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A040A

SPINE HANDPIECE, L300, STAINLESS STEEL **BALLTIP BURR IN SUPPORT TUBE**

Ball diameter: ø 4mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A050A

SPINE HANDPIECE, L220, STAINLESS STEEL **BALLTIP BURR IN SUPPORT TUBE**

Outer diameter: Ø 4mm

Ball diameter: ø 5mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A060A

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 5mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A070A

SPINE HANDPIECE, L220, STAINLESS STEEL **BALLTIP BURR IN SUPPORT TUBE**

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A080A

SPINE HANDPIECE, L300, STAINLESS STEEL **BALLTIP BURR IN SUPPORT TUBE**

Outer diameter: ø 4mm

Ball diameter: ø 6mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G040A

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 3mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G060A

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 3mm x 5mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G040B

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G060B

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G040C

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Sterile, Single use (use directly with NSS-120 Motor)

Tip size: ø 5mm x 5mm

NLB-Z-G060C

SPINE HANDPIECE, L300, STAINLESS STEEL

CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 5mm x 5mm

Working length: 300mm Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-D040A

SPINE HANDPIECE, L220, STAINLESS STEEL EGG BURR IN SUPPORT TUBE

Tip size: ø 3mm x 5mm

Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-D060A

SPINE HANDPIECE, L300, STAINLESS STEEL EGG BURR IN SUPPORT TUBE

Tip size: ø 3mm x 5mm

Working length: 300mm Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-D040B

SPINE HANDPIECE, L220, STAINLESS STEEL EGG BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Working length: 220mm Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-D060B

SPINE HANDPIECE, L300, STAINLESS STEEL

EGG BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Working length: 300mm Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-D040C

SPINE HANDPIECE, L220, STAINLESS STEEL

EGG BURR IN SUPPORT TUBE

Tip size: ø 5mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-D060C

SPINE HANDPIECE, L300, STAINLESS STEEL EGG BURR IN SUPPORT TUBE

Tip size: ø 5mm x 5mm Working length: 300mm







NLB-Z-A210B

SPINE HANDPIECE, L220, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A220B

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm Working length: 300mm Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A010B

SPINE HANDPIECE, L220, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A020B

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G140B

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Tip size: ø 3mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G160B

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Tip size: ø 3mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A210C

SPINE HANDPIECE, L220, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm Working length: 220mn Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A220C

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Outer diameter: ø 4mm Ball diameter: ø 3mm Working length: 300mm Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G140C

SPINE HANDPIECE, L220, STAINLESS STEEL EGG BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Tip size: ø 3mm x 5mm Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G160C

SPINE HANDPIECE, L300, STAINLESS STEEL EGG BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Tip size: ø 3mm x 5mm Working length: 300mn

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A010C

SPINE HANDPIECE, L220, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A020C

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G140C

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Tip size: ø 3mm x 5mm Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-G160C


SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH BEVELED PROTECTION



Tip size: ø 3mm x 5mm



STERILE R



For Interlaminar APPROACH SURGERY,

WITH THE HANDPIECES FOR STENOSIS SURGERY, WHILE

SURGERIES.

THE ENDOSCOPE WITH THE LONGER WORKING LENGTH AND

A SHORTER AND A LONGER WORKING LENGTH ENDOSCOPE CAN BE SELECTED TO ENSURE COVERAGE OF THE ENTIRE SPINAL COLUMN.

THE ENDOSCOPE WITH THE SHORTER WORKING LENGTH IS COMPATIBLE

THE HANDPIECES ARE ALL SUITABLE FOR TRANSFORAMINAL APPROACH



NLB-Z-A210A

SPINE HANDPIECE, L220, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 3mm Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A230A

SPINE HANDPIECE, L220, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 4mm Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A010A

SPINE HANDPIECE, L220, STAINLESS STEEL BALL-TIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 3mm Working length: 220mm Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G040B

SPINE HANDPIECE, L220, STAINLESS STEEL CYLIN-DRICAL BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A210B

SPINE HANDPIECE, L220, DIAMOND **BALLTIP BURR IN SUPPORT TUBE** WITH DISTAL PROTECTION

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A010B

SPINE HANDPIECE, L220, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Outer diameter: ø 4mm Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G140B

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Tip size: ø 3mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A210C

SPINE HANDPIECE, L220, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm Working length: 220mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A010C

SPINE HANDPIECE, L220, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G140C

SPINE HANDPIECE, L220, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Tip size: ø 3mm x 5mm







INTERLAMINAL



NLB-Z-A220A

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 3mm Working length: 300mn

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A240A

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 4mm

Working length: 300mm Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A020A

SPINE HANDPIECE, L300, STAINLESS STEEL **BALLTIP BURR IN SUPPORT TUBE**

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G060B

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A220B

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A020B

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G160B

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Tip size: ø 3mm x 5mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A220C

SPINE HANDPIECE, L300, DIAMOND **BALLTIP BURR IN SUPPORT TUBE WITH** BEVELED PROTECTION

Ball diameter: ø 3mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A020C

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH **BEVELED PROTECTION**

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G160C

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Outer diameter: ø 4mm Tip size: ø 3mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A220A

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Ball diameter: ø 3mm Working length: 300mn

Sterile, Single use (use directly with NSS-120 Motor)

NLB-Z-A240A

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE

Outer diameter: ø 4mm Ball diameter: ø 4mm Working length: 300mr

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A020A

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP **BURR IN SUPPORT TUBE**

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G060B

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE

Tip size: ø 4mm x 5mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A220B

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Outer diameter: Ø 4mr Ball diameter: ø 3mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A020B

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Ball diameter: ø 3mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G160B

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH DISTAL PROTECTION

Tip size: ø 3mm x 5mm Working length: 300mn

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A220C

SPINE HANDPIECE, L300, DIAMOND BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm Working length: 300mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-A020C

SPINE HANDPIECE, L300, STAINLESS STEEL BALLTIP BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Ball diameter: ø 3mm

Sterile, Single use (use directly with NSS-120 Motor)



NLB-Z-G160C

SPINE HANDPIECE, L300, STAINLESS STEEL CYLINDRICAL BURR IN SUPPORT TUBE WITH BEVELED PROTECTION

Tip size: ø 3mm x 5mm



48 Product Catalogue 49

WE KNOW THAT

Patient Safety

STARTS AT ENGINEERING.

THE SKILL LEVEL OF A HIGH-SPEED MOTOR SYSTEM IS GIVEN BY THE VARIETY OF
ATTACHMENTS AND THEIR EXTENSIVE FUNCTIONALITY. THE MAIN PRINCIPLES OF OUR
CONSTRUCTION ARE SAFETY, ERGONOMICS AND SUSTAINABILITY.



NeuroSmart Perforator Handpiece

EXPERIENCE THE DIFFERENCE A Handpiece Can Make.

WITH FULL STAINLESS STEEL FINISH, THE NEW GENERATION emd® PERFORATOR HANDPIECE TAKES CRANIAL PERFORATION TO THE NEXT LEVEL. A HIGHLY EFFICIENT YET EASY USE MADE POSSIBLE BY THE INTEGRATED SMART TECHNOLOGY.

The speed reducer features Automatic Trepan Mode recognition for the NeuroSmart Control Unit. Ergonomic handle for a comfortable, stable grip. Handpiece is available with Hudson end connection.

Thanks to the One-Click Connection, the Perforator Handpiece can be simply mounted on the motor. The perforator lock/release knob can be easily operated with the thumb.

The ergonomic handle is very comfortable to hold and its surface finish allows a secure grip.







Illustration of **emd**®Perforator Handpiece equipped with NeuroLine Single Use Cranial Perforator.



Illustration of $\textbf{emd}^{\$} \text{Perforator Handpiece equipped with}$ NeuroLine Reusable Cranial Perforator.



NSH-210 PERFORATOR HANDPIECE



SPEED SET AT 1.000 RPM

emd® PERFORATOR HANDPIECE ACCEPTS HUDSON END CONNECTION

Features to be used with emd his motor and cranial perforators

Unique Advantages. Multiple Safety Mechanisms.

emd NEUROLINE CRANIAL PERFORATORS WERE DESIGNED AND DEVELOPED TO PROVIDE MAXIMUM PERFORMANCE WITH MAXIMUM PATIENT SAFETY. THE INTEGRATED THREE-LEVEL SAFETY TECHNOLOGY IS PATENT PROTECTED AND OFFERS EXCLUSIVE DRILLING EXPERIENCE FOR THE USER.

EXCELLENT TIP CONTROL TO AVOID SKIDDING

The start of perforation is a critical point, the guiding tip of the drill ensures keeping the perforator exactly on the operation point to avoid skidding on uneven surfaces as well.



HIGH PERFORMANCE & CONSISTENT DRILLING QUALITY

Thanks to the specially developed edge design, the simultaneous movement of the drilling and cutting heads guarantees optimum results.



RESTARTABLE AT ANYTIME

The surgeon has the possibility to continuously supervise the perforation, as the device can be stopped and safely restarted anytime.



BONE PAD FORMATION

Cranial perforation is not a simple drilling, NeuroLine perforators form a bone pad at the end of the drilling process, which protects the dura.



Product Versatility & Compatibility

PRODUCT RANGE OF **emd**® NEUROLINE CRANIAL PERFORATORS GIVES A FREE CHOICE FOR THE CUSTOMER EITHER THEY PREFER SINGLE USE OR REUSABLE MODELS.

All our perforators are manufactured with universal Hudson end coupling.

Thanks to the standardized connection, compatibility is granted for the most commonly used medical high speed motordrill trepan adapters/speed reducers, which are equipped with standard Hudson chuck fitting.








IMPORTANCE OF CORRECT POSITIONING

To ensure a safe and successful operation, maintaining the correct position (90°) of the perforator against the skull is crucial and must be remembered at all times when drilling is being performed.

Explanation of Article Numbers by Size




NLM-11/14-3.0 NLM-B/A-C

(eq. for disposable perforator for thin bone (eg. for reusable perforator for thick bone

A: DIAMETER OF PERFORATOR EXTERNAL DRILL PART • DIAMETER OF UPPER PART OF HOLE

B: DIAMETER OF PERFORATOR INTERNAL DRILL PART • DIAMETER OF LOWER PART OF HOLE

C: MINIMUM THICKNESS OF SKULL BONE REQUIRED FOR OPERATION OF PERFORATOR (OVERHANG OF INTERNAL DRILL PART ON EXTERNAL DRILL PART)



SAFETY BEING OUR TOP PRIORITY, SINGLE-USE PERFORATORS ARE EQUIPPED WITH THREE SAFETY SYSTEMS.

THE BUILT-IN SAFETY MECHANISM ALLOWS THE SAME PRES-SURE TO BE APPLIED ON SOFTER AREAS AS WELL WITH SAME EFFICIENCY.

















FOR USE ON SKULLS MIN. 3.0MM THICK



NLO-6/9-3.0

CRANIAL PERFORATOR, SINGLE USE

NLO-7/11-3.0

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile

NLO-9/13-3.0

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile

NLO-11/14-3.0

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile



NLO-6/9-1.5

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile

NLO-7/11-1.5

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile



NLO-9/13-1.5

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile



NLO-11/14-1.5

CRANIAL PERFORATOR, SINGLE USE Disposable, sterile

OUR REUSABLE CRANIAL PERFORATORS WITH PATENTED TECHNOLOGY OFFER HIGHEST SAFETY FOR BOTH USER AND PATIENT.

WITH THE DEVELOPMENT OF THE REUSABLE PERFORATOR, WE AIMED TO CREATE SUCH A REUSABLE DEVICE THAT IS ENTIRELY CLEANABLE AND STERILIZABLE.

CONSISTING OF FOUR PARTS, PRODUCT ASSEMBLY AND DISASSEMBLY IS PARTICULARLY EASY. THE **emd**® REUSABLE CRANIAL PERFORATOR IS CURRENTLY THE ONLY CRANIAL PERFORATOR WITH SUCH PRODUCT FEATURES ON THE MARKET.

AVAILABLE WITH HUDSON CONNECTION.














NLM-6/9-3.0

CRANIAL PERFORATOR, REUSABLE Autoclavable available in set with Holder Device



CRANIAL PERFORATOR, REUSABLE Autoclavable available in set with Holder Device

NLM-11/14-3.0

CRANIAL PERFORATOR, REUSABLE Autoclavable available in set with Holder Device



NLM-S2-6/9 HOLDER DEVICE

Autoclavable Use with perforator NLM-6/9-3.0

NLM-S2-7/11 HOLDER DEVICE

Autoclavable Use with perforator: NLM-7/11-3.0

NLM-S2-11/14

HOLDER DEVICE Autoclavable Use with perforator: NLM-11/14-3.0



NeuroLine Reusable Perforators are fully dismountable. To make the assembly and reprocessing easier, we developed a holder device for storage and sterilization.

Handling is made easy with the components' dedicated locations identified by pictograms. All cutting edges are protected, minimizing the risk of damage to inserted parts.





Accessories emd 59

WE KNOW THAT THE

Success Of An Operation

RELIES ON THE TEAM BEHIND.

THE LIFE EXPECTANCY OF A HIGH QUALITY MEDICAL DEVICE DEPENDS ON BOTH THE ACCURACY OF MANUFACTURING AND THE CARE TAKEN IN ITS OPERATION AND REPROCESSING.

AS MANUFACTURERS, WE STRIVE TO OFFER SOLUTIONS FOR MEDICAL PROFESSIONALS TO ENSURE EASE OF USE, STORAGE AND STERILISATION.



TO KEEP THE OPERATION AREA COOL WE RECOMMEND TO CHOOSE A SUFFICIENT IRRIGATION METHOD.

USING THE INTEGRATED PUMP OF THE NEUROSMART CONTROL UNIT WITH THE SINGLE USE IRRIGATION TUBE SET ADAPTED ON THE HANDPIECE PROVIDES A COMFORTABLE SOLUTION.

THE REUSABLE IRRIGATION ADAPTER CAN BE EASILY AND SAFELY FIXED ON THE HANDPIECE. WIDE SELECTION IS AVAILABLE FOR THE ALL WORKING LENGTH RANGE.



NSI-641

STANDARD IRRIGATION ADAPTER, L25 Standard – Size 1 NSH-421

NSI-642

STANDARD IRRIGATION ADAPTER, L45 Standard – Size 2 NSH-412 NSH-422

NSI-643

STANDARD IRRIGATION ADAPTER, L70 Standard – Size 3 NSH-413 NSH-423

NSI-644

STANDARD IRRIGATION ADAPTER, L95 Standard – Size 4 NSH-414 NSH-424



STANDARD IRRIGATION ADAPTER, L120 Standard – Size 5 NSH-425

NSI-646

STANDARD IRRIGATION ADAPTER, L170 Standard – Size 6 NSH-426

NSI-647

STANDARD IRRIGATION ADAPTER, L220 Standard – Size 7 NSH-427



FLUID BOTTLE HOLDER

NSI-610

IRRIGATION SET made by ASSAmed GmbH / REF AMP.85.3000





The Importance Of **Product Care**

PROPER PRE-CLEANING IS ESSENTIAL TO PROVIDE EFFECTIVE STERILIZATION. THANKS TO OUR WASH-RINSE DEVICE ALL THE ATTACHMENTS ARE FULLY FLUSHABLE.

THE SYSTEM COMPONENTS CAN BE EASILY POSITIONED AND ITS USAGE ENSURES THAT THE INSIDE OF EACH COMPONENT IS PERFECTLY FREE OF CONTAMINATION.



fix the system elements safely, granting stable position during the preprocessing.

> We strongly recommend using a synthetic Oiling Spray dedicated for cleaning and lubrication, as well as following our product specific

> > instructions of use.

NSC-710

WASH-RINSE DEVICE with cleaning connections for attachments Autoclavable



NSC-712

CERAMIC FILTER For NSC-710




NSC-720

OILING SPRAY Use for all handpieces

NSC-721

OILING SPRAY NOZZLE



NSC-730

INSTRUMENT TRAY with lid and silicon holders Autoclavable dimensions: 482mm x 246mm x 73mm (L x W x H) For storage and sterilization

Notes		Notes	

Notes	



emd-med.com

EMD KFT. H-4031 Debrecen, Bartók Béla u. 113/B +36 52 486 034, +36 52 422 465 info@emd.hu

2023_EN_V2